SimplifyMedical

The Difference ISClear®



Physiologic- Variable Center of Rotation

Radiologic- MRI Compatible*





The Difference is Clear

ANATOMIC: Lowest Available Disc Heights, Starting at 4mm

Designed to fit patient anatomy, protect facets and prevent overstuffing.



4mm




5_mm



5mm


PHYSIOLOGIC: Variable Center of Rotation


3-piece semi-constrained design to replicate natural motion of a healthy disc.



RADIOLOGIC: MRI Compatible*

Detailed anatomic visualization, minimal artifact and radiation-sparing.









Simplify_sMedical

Simplify Medical, Inc. 685 N Pastoria Avenue Sunnyvale, CA 94085 Tel: 650-947-3472 ext. 3 info@simplifymedical.com www.simplifymedical.com

DISCLAIMER: This content is intended for healthcare professionals. For full safety information please visit: www.ifu.simplifymedical.com to view Simplify Disc Instructions for Use, Summary of Safety and Effectiveness and Surgical Technique Guide.

Patented / Patents Pending: www.simplifymedical.com/patents © 2020 Simplify Medical, Inc. Simplify Medical, Simplify Disc and The Difference is Clear are trademarks of Simplify Medical.